
STANDARD Q COVID-19 Ag Saliva


STANDARD™ Q COVID-19 Ag Saliva Test
 PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST SD BIOSENSOR

KIT CONTENTS


SPECIMEN COLLECTION

- **Methods of specimen collection - Saliva with mucus**
- Do not eat, drink, smoke or chew gum for 30 minutes before giving your saliva specimen.
 - The specimen collection cup should not be touched the mouth.
 - The saliva with snorted nasal mucus should be collected enough to cover the bottom of the specimen collection cup.


- Instruct patient to drain mucus from the back of the nose and throat at least 3 times.
- Ask patient to lift specimen collection cup close to his/her mouth and take a deep breath in and cough hard directly into specimen collection cup.
- Instruct patient to lift specimen collection cup close to his/her mouth and spit out the saliva with snorted nasal mucus into specimen collection cup.

PREPARATION AND TEST PROCEDURE

- **Preparation**
- Allow the test device and collected specimen to room temperature prior to testing.
 - Carefully read the instruction for using the STANDARD Q COVID-19 Ag Saliva Test.
*Check the expiry date at the back of the foil pouch. Do not use the kit, if expiry date has passed.
 - Check the test device and the desiccant pack in the foil pouch.

Test Procedure

- Insert the soaked swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
- Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
- Assemble the nozzle cap tightly with the tube.

- Apply 4 drops of extracted specimen to the specimen well of the test device.
 - Read the test result in 15-30 minutes.
-
- Do not read test results after 30 minutes. It may give false results.

CONTROL PREPARATION AND TEST PROCEDURE

- **Positive/Negative control**
- **Preparation**
- Put the STANDARD COVID-19 Ag Positive or Negative Control swab into an extraction buffer tube. While squeezing the buffer tube, stir the swab more than 5 times.
 - Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.
 - Press the nozzle cap tightly onto the tube.
 - Apply 3 drops of extracted specimen to the specimen well of the test device.
-

- Read the test result in 15-30 minutes.
-
- Do not read test results after 30 minutes. It may give false results.

INTERPRETATION OF TEST RESULT

■ **Visual inspection**

Test result	Example	Description
Negative		1. A colored band will appear in the top section of the result window to show that the test is working properly. This band is control line (C).
Positive		2. A colored band will appear in the lower section of the result window. This band is test line of SARS-CoV-2 antigen (T). 3. Even if the control line is faint, or the test line isn't uniform, the test should be considered to be performed properly and the test result should be interpreted as a positive result.
Invalid		

* The presence of any line no matter how faint the result is considered positive.
 * Positive results should be considered in conjunction with the clinical history and other data available.

EXPLANATION AND SUMMARY

■ **Introduction**
 Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "SARS-CoV-2 (COVID-19)", was discovered because of Wuhan Viral Pneumonia cases in 2019, and was named by the World Health Organization on January 12, 2020, confirming that it can cause colds and the Middle East Respiratory Syndrome (MERS) and more serious diseases such as acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.

■ **Intended use**
 STANDARD Q COVID-19 Ag Saliva Test is a rapid chromatographic immunoassay for the qualitative detection of specific antigens to SARS-CoV-2 present in human saliva specimen. This test is for administration by healthcare workers and labs

only, as an aid to early diagnosis of SARS-CoV-2 infection in patient with clinical symptoms with SARS-CoV-2 infection. It provides only an initial screening test result. This product is strictly for medical professional use only and not intended for personal use. The administration of the test and the interpretation of the results should be done by a trained health professional. The result of this test should not be the sole basis for the diagnosis; confirmatory testing is required.

■ **Test principle**
 STANDARD Q COVID-19 Ag Saliva Test has two pre-coated lines, "C" Control line, "T" Test line on the surface of the nitrocellulose membrane. Both the control line and test line in the result window are not visible before applying any specimen. Mouse monoclonal anti-SARS-CoV-2 antibody is coated on the test line region and mouse monoclonal anti-Chicken IgY antibody is coated on the control line region. Mouse monoclonal anti-SARS-CoV-2 antibody conjugated with color particles are used as detectors for SARS-CoV-2 antigen device. During the test, SARS-CoV-2 antigen in the specimen interact with monoclonal anti-SARS-CoV-2 antibody conjugated with color particles making antigen-antibody color particle complex. This complex migrates to the membrane via capillary action until the test line, where it will be captured by the mouse monoclonal anti-SARS-CoV-2 antibody. A colored test line would be visible in the result window if SARS-CoV-2 antigens are present in the specimen. The intensity of colored test line will vary depending upon the amount of SARS-CoV-2 antigen present in the specimen. If SARS-CoV-2 antigens are not present in the specimen, then no color appears in the test line. The control line is used for procedural control, and should always appear if the test procedure is performed properly and the test reagents of the control line are working.

■ **KIT STORAGE AND STABILITY**
 Store the kit at 2-30°C / 36-86°F out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

■ **QUALITY CONTROL**
 A control line is used in the test as a procedural control. A visible control line confirms that the lateral flow of the test is successful but is not the confirmation that the specimen and buffer have been applied properly.

- **External quality control**
- Positive and negative controls may be supplied with each kits or available separately from distributor of SD BIOSENSOR.
 - Positive and negative controls should be performed as a real specimens.
 - It is recommended that positive and negative controls be run:
 - once for each new lot.
 - once for each untrained operator.
 - as required by test instructions and in accordance with local, states and federal regulations of accreditation requirements.

STANDARD COVID-19 Ag Positive Control swab		
Result: T Line (SARS-CoV-2 antigen)	Interpretation	Follow up
Positive	PASS	-
Negative	FAIL	Retest
No Control (C) Line	Invalid	Retest

STANDARD Respiratory Negative Control swab		
Result: T Line (SARS-CoV-2 antigen)	Interpretation	Follow up
Negative	PASS	-
Positive	FAIL	Retest
No Control (C) Line	Invalid	Retest

If either the control invalid result recurs, contact SD BIOSENSOR Technical Support.

WARNINGS AND PRECAUTIONS

- Bring the kit contents and the specimens to room temperature before testing.
- Do not re-use the test kit.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not use the extraction buffer tube of another lot.
- Do not smoke, drink or eat while handling specimen.
- Wear personal protective equipment, such as gloves and lab coats when handling kit reagents. Wash hands thoroughly after the tests are done.
- Clean up spills thoroughly using an appropriate disinfectant.
- Handle all specimens as if they contain infectious agents.
- Observe established precautions against microbiological hazards throughout testing procedures.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste. Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.
- Desiccant in foil pouch is to absorb moisture and keep humidity from affecting products. If the moisture indicating desiccant beads change from yellow to green, the test device in the pouch should be discarded.

LIMITATION OF TEST

- The test procedure, precautions and interpretation of results for this test must be followed strictly when testing.
- The test should be used for the detection of SARS-CoV-2 antigen in human saliva specimens only, other specimen types have not been validated.
- This test can not be used for quantifying SARS-CoV-2 antigen concentration.
- Failure to follow the test procedure and interpretation of test results may adversely affect test performance and/or produce invalid results.
- The test result must always be evaluated with other data available to the physician.
- A negative result may occur if the concentration of antigen in a specimen is below the detection limit of the test or if the specimen was collected or transported improperly, therefore a negative test result does not eliminate the possibility of SARS-CoV-2 infection, and should be confirmed by viral culture or molecular assay.
- Positive test results do not rule out co-infections with other pathogens been validated with the assay.

[Malaysia result]
 Table 2. STANDARD Q COVID-19 Ag Saliva Test performance. RT-PCT test was used as the comparator method in the study.

STANDARD Q COVID-19 Ag Saliva Test	PCR		
	Positive	Negative	Total
	Positive	18	0
Negative	1	73	74
Total	19	73	92
Sensitivity (N, 95% CI)	94.74% (18/19, 73.97% - 99.87%)		
Specificity (N, 95% CI)	100% (73/73, 95.07% - 100%)		

Table 3. STANDARD Q COVID-19 Ag Saliva Test performance against the PCR. - by Cycle Threshold

STANDARD Q COVID-19 Ag Saliva Test	PCR	
	Positive (Ct ≤ 33)	Positive (Ct > 33)
	Positive	16
Negative	0	1
Total	16	3
Sensitivity (95% CI)	100% (79.41% - 100%)	33.33% (0.84% - 90.57%)

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PERFORMANCE CHARACTERISTICS

■ Limit of Detection (LoD)

The SARS-CoV-2 positive specimen was prepared by spiking inactivated SARS-CoV-2 (2019-nCoV) NCCP 43326/2020/Korea strain to SARS-CoV-2 negative saliva specimen confirmed with PCR. LoD is determined as 3.47 X 10⁻² TCID₅₀/ml.

■ Clinical evaluation

The prospective diagnostic evaluation of STANDARD Q COVID-19 Ag Saliva Test with a total number of enrolled individuals of 101 was conducted in Brazil and Malaysia.


[Brazil result]
 Table 1. STANDARD Q COVID-19 Ag Saliva Test performance within 8 days of symptom onset specimen. RT-PCR test was used as the comparator method in the study.

STANDARD Q COVID-19 Ag Saliva Test		PCR		
		Positive	Negative	Total
		5	0	5
Negative	0	4	4	
Total	5	4	9	
Sensitivity (N, 95% CI)	100% (5/5, 47.82% - 100%)			
Specificity (N, 95% CI)	100% (4/4, 91.59% - 100%)			

STANDARD Q COVID-19 Ag Saliva
STANDARD™ Q COVID-19 Ag Saliva Test

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

BIBLIOGRAPHY

- Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) infection is suspected. Interim guidance. WHO.2020
- Diagnostic detection of Wuhan coronavirus 2019 by real-time RT-PCR.2020
- Diagnosis and treatment of pneumonia caused by new coronavirus (trial version 4) National Health Commission. 2020
- Guo L et al. Profiling Early Humoral Response to Diagnose Novel Coronavirus Disease (COVID-19). Clinical Infectious Disease. 2020
- Zhao J et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. Clinical Infectious Disease. 2020
- How to safely collect sputum samples from patients suspected to be infected with pneumonic plague. WHO Interim guidance. September 2016.

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SYMBOL

	Reference number		Caution		Use by		Batch code		Consult Instructions for Use		Do not re-use
	In vitro Diagnostics		Note		Manufacturer		Date of manufacture		Contains Sufficient for <=9 Tests		Keep away from sunlight
	Indicate that you should keep the product dry		To indicate the temperature limitations in which the transport package has to be kept and handled.		CE		Fullfill the requirements of Directive 98/79/EC on in vitro diagnostic medical devices		Do not use if packaging is damaged		

Authorized Representative
MT Promed Consulting GmbH
Altenhofstrasse 80 66386 St. Ingbert Germany
Phone : +49 6894 581020, Fax : +49 6894 581021

Manufactured by **SD Biosensor, Inc.**
Head office : C-4th&5th, 16, Deogyong-daero 1556beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, 16690, REPUBLIC OF KOREA
Manufacturing site : 74, Osongsaemyeong 4-ro, Osong-eup, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, 28161, REPUBLIC OF KOREA

Please contact us for any complaints/inquiries/suggestions via email (sales@sdbiosensor.com), phone (+82-31-300-0400) or website (www.sdbiosensor.com).